

POLSTERUNG

UPHOLSTERY

Die Polsterung der Untersuchungsliegen, Hocker und Stühle bestehen aus einer mit Schaumstoff verklebten, stabilen Trägerplatte. Diese ist nahtlos mit einem strapazierfähigen, reißfesten und für den medizinischen Einsatz geeigneten Kunstlederbezug versehen. Es steht ein reichhaltiges Spektrum an modernen Polsterfarben zur Auswahl.

Abgerundete und verstärkte Polsterkanten vervollständigen das edle und funktionelle Design.

Die Kopfteilverstellung und -fixierung der Untersuchungsliegen erfolgt mittels stabiler, korrosionsgeschützter Verstellraster.

The upholstery of the examination couches, stools and chairs consists of a foam padding glued on a synthetic material board. It is seamlessly covered with a tearproof synthetic leather covering which is suitable to medical application. A wide range of modern upholstery colours is available.

Rounded edges and strong cushions complete the noble and functional design.

The head section adjustment and fixation of the couch is carried out by a stable galvanized metal-grid.

Gestellfarben/ frame colours:



Kunststoffbeschichtung:
Powder-coating:
Weiß (RAL 9010)
White (RAL 9010)



Kunststoffbeschichtung:
Powder-coating:
Weißaluminium (RAL 9006)
White-aluminium (RAL 9006)



Verchromt
Chromium-plated

Polsterfarben / upholstery colours:

	-50	Weiss White
	-56	Beige Beige
	-67	Limonade Lemonade
	-75	Agave Agave
	-73	Apricot Apricot
	-53	Perle Pearl
	-58	Stone Stone
	-55	Dolphin Dolphin
	-62	Bleu Light blue
	-74	Atoll Atoll
	-80	Royal Royal
	-107	Graphitgrau Graphite grey
	-83	Schwarz Black
	-96	Schwarz, elekt. leitend Black, electrically conducting